Stretch Film

ExpressWrap (Extended Core Hand Wrap)

Values (listed by film thickness) Test 80 gauge **Reference Methods** 20.3 micron 0.80 mil **Ultimate Tensile MD** (psi) ASTM D882 6300 **Ultimate Tensile TD** (psi) ASTM D882 4250 **Elongation MD** (%) ASTM D882 600 **Elongation TD** (%) ASTM D882 800 **Tear Resistance MD** 260 grams **ASTM D1922** (grams and g/mil) 325 g/mil 500 **Tear Resistance TD** grams **ASTM D1922** (grams and g/mil) 625 g/mil **Dart Drop** 150 grams ASTM D1709 (grams and g/mil) 187.5 g/mil 250 **Pre-stretch** (% at break) IPG/Highlight Clina IPG method 150 grams (in to out, grams at 200% pre-stretch)

While we believe them to be reliable, the statements and information herein are only for general guidance and are not warrants or guarantees for accuracy and completeness. The user must, by test or otherwise, determine suitability for this purpose. There is no warranty of fitness for a particular purpose. Our standard term and conditions of sale apply exclusively to all orders, and all liability for damages of any kind, including consequential, exceeding purchase price is excluded. No one is authorized by us to make oral warranties. We reserve the right to make changes without notice or obligation in our products and publications.



100 Paramount Drive, Suite 300 | Sarasota, FL 34232 | USA Customer Service: 800.IPG.8273 | 800.474.8273

Tape Technical Service: 877.447.4832 www.itape.com | info@itape.com

EFFECTIVE: 04/14



